


Photoactivatable and Membrane permeant Ins(1,4,5)P₃ - caged IP₃

Membrane-permeant derivatives of inositol polyphosphates require the intracellular enzymatic hydrolysis of several protecting groups, for instance of acyloxymethyl esters, in order to generate the biologically active compound. The highly complicated kinetics of these biochemical steps may lead to unphysiological effects. The physiological signal usually appears to be very rapid. The photolysis of membrane-permeant caged derivatives of


Ins(1,4,5)P₃ mimic fast intracellular responses. In an initial step cells are loaded with the caged Ins(1,4,5)P₃/PM derivative. Within 30-180minutes all bioactivatable protecting groups remove, generating caged inositol polyphosphate. The cage is known to prevent biological activity when placed at the right position, in this case the 6-hydroxy-group.



The photochemical destruction of the cage (~360nm) releases active Ins(1,4,5)P₃ within a few seconds, thus mimicking the rapid responses of the receptor / phospholipase C signaling system in the cell (1).


Experiments

Since this approach does not directly trigger other signaling events, for instance PKC activation after receptor-mediated diacylglycerol (DAG) formation, membrane-permeant derivatives of signaling molecules are able to help dissecting signaling pathways. In the shown experiment HeLa cells were co-loaded with 2μM caged Ins(1,4,5)P₃/PM [cag-iso-2-145] and Oregon Green 488 BAPTA-1AM. All compounds were loaded for 80min at room temperature. Stimulus was 1 burst of 5 UV flashes (300-400nm bandwidth) delivered 42s after the onset of the recording: [Ca²⁺]_i raised! (2)



2

Experiments performed by Dr. Valeria Piazza and analysed by Dr. Catalin D. Ciubotaru in the laboratory of Prof. Fabio Mammano at the Venetian Institute of Molecular Medicine – Padua University, Italy, <http://www.vimm.it>



D-2,3-O-Isopropylidene-6-O-(2-nitro-4,5-dimethoxy)benzyl-*myo*-Inositol 1,4,5-trisphosphate-Hexakis(propionoxymethyl) Ester

Product No.: cag-iso-2-145-10 (10*10µg)
cag-iso-2-145-100 (1*100µg)

includes 200µL Pluronic® F-127 in DMSO (10%)

Formula: C₄₂ H₆₄ N O₃₁ P₃
MW: 1171.27

Preparation

It is probably a good idea to aliquot the sample. The compound is soluble in CH₂Cl₂ or DMSO, which evaporates instantly under reduced pressure. The evaporation vessel should be filled with argon (better) or nitrogen afterwards. The compound is sensitive to water on a longer time scale. Therefore, please store the compound in substance or in dry DMSO (for not longer than 2 weeks) at -20°C or below. The freezing process should be performed very quickly (-80°C), not just in the freezer. For incubations, dissolve an aliquot of the cell penetrating compound in dry DMSO.

Take out a small amount (e.g. 1µl) and mix with same amount of Pluronic® F127 in DMSO (10%).

To this mixture add 100µl of the **serum-free cell supernatant**, mix thoroughly with a pipette and immediately add back to the cells. The final DMSO concentration in the experiment should not exceed 0.5%. The final concentration of caged-iso-Ins(1,4,5)P₃/PM should be in the **1-3µM** range, depending on the cell type. For calcium measurements after de-caging use one of the standard calcium sensors. **Incubation** in the dark at room temperature for the cell penetrating compound and the calcium sensor should be **30-120min**. Subsequently, it is possible to return to different buffers (with serum, if necessary). If buffer is not changed, **calcium levels** can be measured **within 5-10min**. If buffers are changed, a longer adaptation phase (30 min) is recommended. To **un-cage** caged-iso-Ins(1,4,5)P₃/PM, scan cells once with an excitation around **360nm** of an argon-ion UV laser or another UV light source

Summary:

- stock solution in CH₂Cl₂ or DMSO (storage: 2 weeks, -20°C)
- final concentration of caged-iso-InsP₃/PM: 1-3 µM
- incubation: 30-120min.
- un-cage caged-iso-InsP₃/PM: with an argon-ion UV laser using the 345-355nm line

Other caged Inositol Phosphates

[cag-iso-2-145]	[cag-6-145]	[cag-0-145]
caged-InsP ₃ -DMNB Membrane-permeant and photolabile derivate of Ins(1,4,5)P ₃ with DMNB (D-23-O-Isopropylidene-6-O-(2-nitro-4,5-dimethoxy)	caged-InsP ₃ -DMNB Photolabile derivate of Ins(1,4,5)P ₃ with the same caged group as the membrane-permeant derivate [cag-iso-2-145]. It is photolyzed with UV light about three times more efficiently than the widely used D- <i>myo</i> -Inositol Trisphosphate-NPE (P ⁴ -1-(2-Nitrophenyl)ethyl Ester) [cag-0-145]	caged-InsP ₃ -NPE Photolabile derivate of Ins(1,4,5)P ₃ with NPE (P ⁴ -1-(2-Nitrophenyl)ethyl Ester)
10 * 10µg / 1 * 100µg	1 * 100µg	1 * 100µg

References

- Adam Bartok et al. ; **NATURE COMMUNICATIONS** 2019. doi.org/10.1038/s41467-019-11646-3
- Lock, J. et al.; **Cell Calcium** 2017 May; 63: 43-47. doi:10.1016/j.ceca.2016.11.006
- Emanuel Oekeke et al., **Biochem. J.** (2016) 473, 757–767 doi:10.1042/BJ20150364
- Vyacheslav M. Shkryl, Joshua T. Maxwell and Lothar A. Blatter: A novel method for spatially complex diffraction-limited photoactivation and photobleaching in living cells
J Physiol 590.5 (2012) pp 1093–1100
- Ian Parker and Ian F. Smith: Recording single-channel activity of inositol trisphosphate receptors in intact cells with a microscope, not a patch clamp
J. Gen. Physiol. 2010, Vol. 136 No. 2 119–127
- Steven M. Wiltgen, Ian F. Smith, and Ian Parker: Superresolution Localization of Single Functional IP3R Channels Utilizing Ca²⁺ Flux as a Readout
Biophysical Journal Volume 99 July 2010 437–446
- Ian F. Smith, Ian Parker: Imaging the quantal substructure of single IP₃R channel activity during Ca²⁺ puffs in intact mammalian cells
PNAS _ April 14, 2009 _ vol. 106 _ no. 15
- Ian F. Smith, Steven M. Wiltgen, Ian Parker: Localization of puff sites adjacent to the plasma membrane: Functional and spatial characterization of Ca²⁺ signaling in SH-SY5Y cells utilizing membrane-permeant caged IP₃;
Cell Calcium 45, 65–76 (2009)
- Hua, X. et al.: Ca²⁺-Dependent Glutamate Release Involves Two Classes of Endoplasmic Reticulum Ca²⁺ Stores in Astrocytes.
J. Neurosci. Res. 76, 86 (2004).
- Samways, D.S.K. et al.: Co-incident signalling between μ -opioid and M3 muscarinic receptors at the level of Ca²⁺ release from intracellular stores: lack of evidence for Ins(1,4,5) P3 receptor sensitisation; **Biochem. J.** 375: 713-720 (2003)
- Wagner, L.I. et al.: Phosphorylation of Type-1 Inositol 1,4,5-Trisphosphate Receptors by Cyclic Nucleotide-dependent Protein Kinases,
J. Biol. Chem. 278 (46): 45811-45817 (2003)
- Li, W.-H. et al.: Cell-permeant caged InsP3 ester shows that Ca²⁺ spike frequency can optimize gene expression;
Nature 392: 936-941 (1998)

For **purchase orders** please go to our online shop or send us an e-mail (info@sichem.de) or fax (+49-421-2208-226)
For any further questions, don't hesitate contacting us !